User Manual (English)

Ver.1.0

Fingertip Pulse Oximeter

This Fingertip Pulse Oximeter is a kind of innovated medical device with non-invasive features for artery SpO2 and PR detection. Being portable, it is able to measure SpO2 and PR values quickly and precisely.

General Description

Haemoglobin Saturation is the percentage between the capacity of Oxihemoglobin (HbO2) that compounded with oxygen and that of all combinative haemoglobin (Hb) in blood. In other words, it is the saturation of Oxihemoglobin in blood. It is a very important physiological parameter for Respiratory and Circulation Systems. Many respiratory diseases could reduce haemoglobin saturation in human blood. Moreover, factors such as Automatic Organic Regulation Malfunction caused by anaesthesia, trauma resulted from major operation and some medical examination can also cause problems in oxygen supply, which might reduce human haemoglobin saturation. As a result, such symptoms as megrim, vomiting and asthenia might appear to patients. Hence, it is very important to know hemoglobin saturation of patient timely in clinical medical aspects.

The fingertip pulse oximeter features in small volume, low power consumption, convenient operation and portability. It is only necessary for patient to put one finger into fingertip photoelectric sensor for diagnosis, and the display screen will directly show measured value of hemoglobin saturation. It has been proved in clinical experiments that it possesses rather high precision and repeatability.

Measurement principle

The principle of the oximeter is as follows: An experience formula of data process is established by exerting Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive hemoglobin and oxihemoglobin (HbO2) in glow and near-infrared zones. Operation principle of the instrument is to combine Photoelectric, Opt., and Infrared Inspection Technology with Capacity Pulse Scanning and Recording Technology, so that two lights with different wavelength (660nm glow and 940nm near infrared light) can be focused on human nail through perspective clamp finger-type sensor. Then measured signal can be obtained by a photosensitive element, information acquired through which will be shown on two groups of LEDs through process in electronic circuits and microprocessor.

11. Operation environment: temperature: 5 C ~ 40 C Humidity: 15% ~ 80%

Atmospheric pressure: 70kPa ~ 106kPa Atmospheric pressure: 70kPa ~ 106kPa

12. Declaration: EMC of this product comply with IEC60601-1-2 standard.

Classification

- 1.Management Class for Medical Devices: I equipment
- 2.Anti-electric Shock Type: Internally powered equipment
- 3.Anti-electric Shock Degree: Type BF equipment

Maintenance and Preservation

- 1.Replace the batteries timely when low voltage lamp is on.
- 2.Clean the surface of fingertip pulse oximeter before it is used to diagnose patients.
- 3.Remove the batteries inside if you will not operate the Oximeter for a long time.
- 4.It would be better to preserve the product in -10~40C (14-104 F) and humidity is 10%~95%.
- 5.It is recommended that the product should be kept dry anytime. A wet ambience might affect its lifetime and even damage the product.
- 6.Please follow the law of the local government to deal with used batteries.

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
Diagram of Operation Principle

- 1. Infrared-ray receiving tube
- 2. Infrared-ray transmitting tube

Precautions for use

- 1.Do not use the Fingertip Pulse Oximeter together with MRI or CT equipment.
- 2.Explosion hazard: Do not use the Fingertip Pulse Oximeter in an explosive atmosphere.
- 3.The Fingertip Pulse Oximeter is intended only as an adjunct in patient assessment. Doctors should make diagnosis in conjunction with clinical manifestation and symptoms.
- 4.Check the Fingertip Pulse Oximeter sensor application site frequently to make sure that the circulation and skin integrity of patient are under good condition.
- 5.Do not stretch the adhesive tape while applying the Fingertip Pulse Oximeter sensor. This may cause inaccurate reading or skin blisters.
- 6.Please read the manual carefully before your operation.
- 7.The Fingertip Pulse Oximeter has no SpO2 prompt, it is not for continuous monitoring.
- 8.Prolonged use or the patient's condition may require changing the sensor site periodically. Change sensor site and check skin integrity, circulatory status, and correct alignment at least every 2 hours.
- 9.Inaccurate measurements may be caused by autoclaving, ethylene oxide sterilizing, or immersing the sensors in liquid.
- 10.Significant levels of dysfunctional hemoglobins (such as carboxyl-hemoglobin or methemoglobin) may cause inaccurate reading.
- 11.Intravascular dyes such as Indocyanine green or methylene blue may cause inaccurate reading.
- 12.SpO2 measurements may be adversely affected in the presence of high ambient light. Please shield the sensor area (with a surgical towel or direct sunlight, for example) if it is necessary.
- 13.Unexpected action may occur if the patient is wearing a defibrillator.
- 14.Medical signal with high frequency or interference caused by defibrillator may lead to inaccurate reading.
- 15.Venous pulsations may cause inaccurate reading.

Product Accessories

- 1.One hang lace
- 2.One user manual

Guidance and manufacture's declaration-electromagnetic radiation-for other EQUIPMENTS and SYSTEMS

1	Guidance and manufacturer's declaration – electromagnetic emission	
2	The Fingertip Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of Fingertip Pulse Oximeter should assure that it is used in such an environment.	
3	Emissions test	Compliance
4	RF emissions CISPR 11	Group 1
5	RF emissions CISPR 11	Class B
6	Harmonic emissions IEC 61000-3-2	N/A
7	Voltage fluctuations/flicker emissions IEC 61000-3-3	N/A

Possible Problems and Resolutions

Problem	Possible reason	Solution
SpO2 or PR can not be shown normally	1. Finger is not plugged correctly 2.Patient's hemoglobin value is too low to be measured	1. Retry by plugging the finger correctly 2.If there is no problem in the product, please go to hospital timely for exact diagnosis
SpO2 or PR is shown unstable	1. The finger might not be plugged deep enough 2. Finger is trembling or the patient is in movement status	1. Retry by plugging the finger correctly 2. Please remain at rest
The Oximeter can not be turned on	1. Inadequate power or power off 2. Batteries might be installed incorrectly 3. The Oximeter might be damaged	1. Please replace the batteries 2. Please reinstall the batteries 3. Please contact with local customer service center
Indication lamps are suddenly off	1. The product automatically shuts off when no signal is detected in 8 seconds 2. Inadequate power	1. Normal 2. Replace the batteries

- 16.It may cause inaccurate reading when the positions of sensor and blood pressure cuff are on the same arterial catheter or intravascular line.
- 17. Hypotension, severe vasoconstriction, severe anemia, or hypothermia may cause inaccurate reading.
- 18.It may cause inaccurate reading by giving use of cardiotoxic to patient after his cardiac arrest or when he is in quiver.
- 19.Bright nail or painted nail may cause inaccurate SpO2 reading. Follow local ordinances and recycling instructions regarding to disposal or recycling of the device and device components, including batteries.

Scope of application / Intended use
The fingertip pulse oximeter can be used to measure human haemoglobin saturation and pulse rate through finger , it can be used in hospitals, families, schools and medical centers.

contraindication: not found

note

- 1.The image in the instruction may have slight differences with the actual instruments.
- 2.The actual parameters and appearance change, without prior notice.

Product include: main machine and SpO2 sensor.
Features
◆ OLED display
◆ Product adopts double color OLED display , can show the six different display mode

- ◆ If the hand movements, under the effect of accelerometer, the interface can save four different kinds of display mode (suitable for matching accelerometer function instrument)
- ◆ Low-power consumption, continuously work for more than six hours with two AAA batteries
- ◆ Low voltage indicator
- ◆ In the absence of signals, the product will be in after 8 seconds to enter standby state
- ◆ Small in volume, light in weight, and convenient to carry
- ◆ Instrument has 5s automatic signal detection function, when you insert finger, timely will automatically start;Automatic startup function instrument (applicable to Automatic startup function instrument)

Operation Instructions

- 1.Install two AAA batteries into battery cassette before covering its cover.
- 2.Plug one finger into rubber hole of the Oximeter (It is best to plug the finger thoroughly) before releasing the clamp with the nail upwards.

Reserves the right to technical change appearance, our products are subject to change without prior notice, please forgive me!

Statement:

- 1, Maintenance with data such as circuit diagram, components list, figure and the detailed rules of correction injection, available only to the repair factory training qualified personnel and units.
- 2, The company can be in the form of email or other electronic files provide users with random files.
- 3, The instrument is not used for evaluation of blood oxygen probe pulse and pulse blood oxygen monitor accuracy.

After-sales service

Ensure that users

- ◆ Please read user manual before using the instrument;
- ◆ According to the requirement of the instruction manual for the operation and daily maintenance, and make sure the machine power supply, and environmental requirements

Maintenance time

Monday to Friday 9 to 17:30pm, except for the national legal holiday

Maintenance regulations

- ◆ To conform to the regulations, free maintenance within the scope of products, with warranty card for free maintenance.All that is beyond the scope of maintenance, provide paid services.
- ◆ With warranty card and shopping invoice,main machine for a year, accessories for three months are under free maintenance services from the date of purchase.
- ◆ Following does not belong to the scope of free maintenance
 - ◆ The fault caused by human factors, the damage;
 - ◆ Due to the use to be inconsistent with the provisions of our company work environment to cause damage;
 - ◆ Due to the product in the our company authorized personnel disassembling or repairing damage;
 - ◆ Products beyond the warranty period;
- ◆ If any problem , please call us in: 9: 00 am to 5: 30 pm from Monday to Friday(except national holiday),call us: 400-828-6667

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- 3.Press button on the front panel; (Note: if equipped with automatic startup function refers to clamp oximeter; need not press the button, the instrument has 5 s automatic signal detection function, directly inserted into the finger, instrument is automatically switched on timely)
- 4.Don't tremble your finger when the Oximeter is working. Your body is not recommended on moving status.
- 5.Press the button on the front panel, if we want change display direction;
- (Note: if equipped with accelerometer function of instrument then don't press the button, hand movements, the instrument with the accelerometer has four corresponding interface switch)
- 6.Read relevant datum from display screen.
- 7. The instrument has the function of sleep, no signal will enter standby state of sleep;
- 8.Please replace new batteries when OLED indicates the batteries are in low power.

When plugging your finger into the Oximeter, your nail surface must be upward.

Declaration: Please use the medical alcohol to clean the rubber before each test and clean the tested finger with alcohol before and after the test. (The rubber inside of the Oximeter adopts medical rubber, which has no toxin, no harm, and brings no side effect such as allergy to the our skin)

Battery installation

- 1. According to the positive and negative identity right will be two AAA batteries into the battery
- 2. In the direction of the arrow at the bottom of the level of battery back cover.

Please note: pay attention to the battery positive and negative polarity, must be installed correctly, otherwise it may cause damage to instrument.

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