	Keep dry		Keep away from sunlight
	Store between 2°C-30°C (36°F-86°F)		Manufacturer
	Expiry date		Date of manufacture
	Recyclable		Batch code
	Catalog number		Do not re-use
	Consult instructions for use		

Distributed by
Easy Healthcare Corporation,
360 Shore Dr., Burr Ridge, IL USA 60527

Email: service@healthcare-manager.com

To learn more, please visit us at www.healthcare-manager.com

Made in China

Check 10 PARAMETERS user manual in more languages

Reagent Strips for Animal Urinalysis

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REF 10 PARAMETERS

10 Parameters Urine Reagent Strips for Animals

Product Name: 10 Parameters Reagent Strips for Animal Urinalysis
Product Description: Easy@Home® 10 Parameters Reagent Strips for Animal Urinalysis are dip-and-read test strips for *in vitro* diagnostic to help detect the presence of Leucocytes, Nitrite, Urobilinogen, Protein, pH, Blood, Specific Gravity, Ketone, Bilirubin, and Glucose in urine. Test results may provide information regarding the status of carbohydrate metabolism, kidney and liver function, acid-base balance, and urinary tract infection. It is measured by comparing the colored reactants attached to a plastic strip with the color chart blocks printed on the bottle label.

Intended Use

Easy@Home® 10 Parameters Reagent Strips for Animal Urinalysis are used for qualitative and semi-quantitative detection of Leukocytes, Nitrite, Urobilinogen, Protein, pH, Blood, Specific Gravity, Ketone, Bilirubin, and Glucose presence in the urine of cats and dogs. Urine undergoes many changes during state of disease or body dysfunction before blood composition is altered to a significant extent. Testing urine is used as an indicator of health or disease, and as such, is a part of routine health screening. Easy@Home® 10 Parameters Reagent Strips for Animal Urinalysis can be used in general evaluation of health, and aids in the diagnosis and monitoring of metabolic or systemic diseases that affect kidney function, endocrine disorders, and diseases or disorders of the urinary tract. For *in vitro* diagnostic use in animals only.

Principle

Leukocytes: Normal urine of dogs and cats specimens generally yield negative results. Trace result suggests a possible kidney infection. When trace results occur, it is recommended to retest using a fresh specimen from the same cat or dog. Repeated trace and positive results are of clinical significance. Always contact your physician with questions or concerns.

However, leukocyte field is often falsely positive in cats. False positive may also occur in specimens infected with E. coli.

Nitrite: Nitrite is not detectable in normal urine of dogs and cats. The presence of infections caused by nitrate-reducing bacteria suggests a possible urinary tract infection. And the nitrite area will be positive. However, the nitrite area will be positive in some cases of infection, depending on how long the urine specimens were retained in the bladder prior to collection. Retrieval of positive cases with the nitrite test ranges from as low as 40% in cases where little bladder incubation occurred, to as high as approximately 80% in cases where bladder incubation took place for at least 4 hours.

Urobilinogen: A small amount of urinary urobilinogen is normal. Urobilinogen, formed from bilirubin by intestinal microflora, is absorbed into the portal circulation and excreted renally. Urobilinogen is one of the major compounds produced in heme synthesis and is a normal substance in urine. Increased urinary urobilinogen occurs with hyperbilirubinemia. The expected range for normal urine with this test is 3.2-16 µmol/L.

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pH: Used to help measure acidity or alkalinity in the body through urine. Normal dogs and cats may have a urine pH of 5.5 to 7.5.

Blood: The tracing result is 10 $\text{caecili}/\mu\text{l}$, and it is required to make clinical judgments based on the specific situation of the specimen. More than 10 $\text{caecili}/\mu\text{l}$ result is abnormal, and a further investigation is required. Often but not necessarily, blood is found in the urine of UTI, kidney infection.

Specific Gravity: Randomly collected urine of dogs may vary in specific gravity from 1.015-1.045. Randomly collected urine of cats may vary in specific gravity from 1.035-1.060. When the urine specific gravity is continuously lower than 1.008, the kidneys may lack the ability to concentrate urine, please judge based on clinical symptoms.

Ketone: Ketones are normally not present in urine of dogs and cats. Ketones are a byproduct of fat break down. When enough of them built up in the blood, it becomes poisonous to the animal. This type of fat break down usually occurs when the animal can't use glucose as a form of energy due to a lack of insulin.

Bilirubin: Most healthy dogs and cats have no bilirubin in their urine. However, because of the low kidney threshold of dogs, a small number of normal dogs contain detectable bilirubin in their urine. Bilirubin can be detected before serum bilirubin is elevated. There is no bilirubin in normal cat urine. Abnormal amounts of bilirubin in the urine are associated with liver disease. False-negative reactions may occur with large amounts of urinary ascorbic acid.

Glucose: When dogs and cats are highly excited, eat too much glucose or fructose, eat a lot of carbohydrate-rich feed, or have kidney disease, such as diabetes, glucose may appear in the urine. In addition, cats may develop temporary hyperglycemia and diabetes after severe stress.

Caution

Water cannot be used as negative quality control liquid. Antiseptic properties of urine cannot prevent the ketone, bilirubin, and urobilinogen from deteriorating. For stale urine specimen, the test results of glucose, pH, nitrite, and blood can be affected because of bacterial growth.

- Do not use after the expiration date.
- Do not touch the reagent areas of the strip.
- Discard any colored strips that may have deteriorated.
- The strips should remain in the closed bottle until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- Used strips should be discarded after testing.

Storage and Stability

Store at 2°C-30°C (36°F-86°F) in the sealed canister up to the expiration date.

Keep out of direct sunlight. The strip is stable through the expiration date printed on the canister label. **DO NOT REMOVE THE DESICCANT.** Remove only enough strips for immediate use. Replace cap immediately and tightly. **DO NOT FREEZE.** Do not use beyond the expiration date.

Note: Once the canister has been opened, the remaining strips are stable for up to 3 months.

Stability may be reduced in high humidity conditions.

Specimen Collection and Preparation

Use only clean dry container to collect urine. Shake the collected urine before testing and test within 2 hours of being collected. Any operations must be in a sanitary environment. The use of urine preservatives is not recommended. Prolonged storage of unpreserved urine at room temperature may result in microbial proliferation.

Materials Provided

- Strips
- Package insert

Materials Needed But Not Provided

Container to collect urine
Timer or watch

Test Conditions

Room temperature: 20°C-30°C (68°F-86°F)

Relative humidity: ≤80%

Test Procedure

1. Remove the strip from the closed canister and use it as soon as possible. Immediately close the canister tightly after removing the required number of strip(s). Completely immerse the reagent areas of the strip in fresh, well-mixed urine and immediately remove the strip to avoid dissolving the reagents. **See illustration 1 below.**

2. While removing the strip from the urine, run the edge of the strip against the rim of the urine container to remove excess urine. Hold the strip in a horizontal position and bring the edge of the strip into contact with an absorbent material (e.g. a paper towel) to avoid mixing chemicals from adjacent reagent areas and/or soiling hands with urine. **See illustration 2 below.**

3. Read the test results carefully within 60-120 seconds in a good light and with the test area held near the appropriate color chart on the bottle label. Changes in color that appear along the edges of the test pads or after more than 2 minutes have passed are of no diagnostic significance.

Results with Leukocytes test portion can be read within 120 seconds.

Nitrate, Urobilinogen, Protein, pH and Blood can be read in 60 seconds.

Specific Gravity in 60 seconds, Ketones in 60-120 seconds and Bilirubin and Glucose in 60 seconds. **See illustration 3 below.**

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