## Multi-Drug Screen Test

- EN -

Apply to the following Cat. No.: EDOAP-754, EDOAP-264

Please refer to the information on the box or pouch for the catalog.

Easy@Home® Multi-Drug Screen Test offers fast, accurate and easy drug testing. Easy@Home® Multi-Drug Screen Test (Cat. No.: EDOAP-754, Cat. No.: EDOAP-264) can detect 5 or 6 of the following commonly abused drugs: Amphetamine (AMP), Benzodiazepines (BZD), Cocaine (COC), Marijuana (THC), Methamphetamine (MET), Opiate (OPI 2000).

This package insert applies to Easy@Home® Multi-Drug Screen Test (Cat. No.: EDOAP-754 and Cat. No.: EDOAP-264). Therefore, some information on the performance characteristics of the product may not be relevant to your test. Please refer to the labels on the packaging and the prints on the test strip to identify which drugs are included in your test.

Easy@Home® Multi-Drug Screen Test is a rapid one step test for the qualitative detection of drug abuse and their principal metabolites in human urine at specified cut off level.

**For in vitro diagnostic use only.**  
**INTENDED USE**

Easy@Home® Multi Drug Screen Test is rapid urine screening test. The test is a lateral flow, one-step immunoassay for the qualitative detection of specific drugs and their metabolites in human urine at the following cut off concentrations:

Drug(Identifier)	Calibrator	Cut-off level	Minimum detection time	Maximum detection time
Amphetamine (AMP)	d-Amphetamine	1000ng/ml	2-7 hours	1-2 days
Benzodiazepine (BZD)	Oxazepam	300 ng/ml	2-7 hours	1-2 days
Cocaine (COC)	Benzoylgegonine	300 ng/ml	1-4 hours	24+ days
Marijuana (THC)	11-nor-d <sub>9</sub> -THC-9-COOH	50 ng/ml	2 hours	Up to 5 days
Methamphetamine (MET)	Methamphetamine	1000ng/ml	2-7 hours	2-4 days
Opiate (OPI2000)	Morphine	2000ng/ml	2 hours	2-3 days

This assay provides only a preliminary test result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary results are positive.

### PRINCIPLE

#### DRUG TESTS

Easy@Home® Multi-Drug Screen Test is a competitive immunoassay that is used to screen for the presence of drugs of abuse in urine. It is chromatographic absorbent device in which drugs or drug metabolites in a sample competitively combine to a limited number of antibody-dye conjugate.

When testing, the urine is absorbed upward by capillary action, mixes with the antibody-dye conjugate, and flows across the pre-coated membrane.

When sample drug levels are at or above the target cutoff, the drug in the sample binds to the antibody-dye conjugate preventing the antibody-dye conjugate from binding to the drug-protein pre-coated in the test region (T). This prevents the development of a distinct colored band in the test region indicating a potentially positive result.

When sample drug levels are zero or below the target cutoff (the detection sensitivity of the test), antibody-dye conjugate binds to the drug-protein pre-coated in the test region (T) of the device. This produces a colored test line that, regardless of its intensity, indicates a negative result.

To serve as a procedure control, a colored line will appear on the control region (C), if the test has been performed properly.

### WARNINGS AND PRECAUTIONS

- This kit is for external use only. Do not swallow.
- Discard after first use. The test cannot be used more than once.
- Do not use test kit beyond expiration date.
- Do not use the kit if the pouch is punctured or not well sealed.
- Keep out of the reach of children.

### STORAGE AND STABILITY

- Store at 4°C-30°C (40°F-86°F) up to the expiration date.
- Keep away from sunlight, moisture and heat.
- DO NOT FREEZE.
- After opening the pouch, please use the drug panel as soon as possible, no more than 1 hour.

### MATERIAL

#### Material provided

- One pouch containing a test panel and a desiccant.
- Package insert

### MATERIAL REQUIRED BUT NOT PROVIDED

- Timer
- Urine cup

### SPECIMEN COLLECTION AND PREPARATION


Collect a urine sample in the urine cup. Urine specimens may be refrigerated 2°C-8°C (36°F-47°F) and stored up to forty-eight hours. For longer storage, freeze the samples at -20°C (-4°F) or below. Bring frozen or refrigerated samples to room temperature before testing. Use only clear aliquots for testing.

### TEST PROCEDURE


Test must be in room temperature 18°C-30°C (65°F-86°F).

1. Open the sealed pouch by tearing along the notch. Remove the test device from the pouch.
2. Hold the one side of the device with one hand. Use the other hand to pull out the cap and expose the absorbent end.
3. Immerse the absorbent end into the urine sample about 10 seconds. Make sure that the urine level is not above the "MAX" line printed on the front of the device.
4. Lay the device flat on a clean, dry, non-absorbent surface.
5. Read the result in 5 minutes. Do not read after 5 minutes.

#### Step 1: Pull the cap off and immerse the strips into urine for 10 seconds


#### Step 2: Read the result in 5 minutes. Do not read after 5 minutes.

### INTERPRETATION OF RESULTS


#### Preliminary Positive (+)

A colored band is visible in each control region(C) . No colored band appears in the appropriate test region(T). It indicates a positive result for the corresponding drug of that specific test region(T).


#### Negative (-)

A colored band is visible in each control region(C) and the appropriate test region(T). It indicates that the concentration of the corresponding drug of that specific test region(T) is zero or below the detection limit of the test.

#### Invalid

If a colored band is not visible in each of the control region(C) or a colored band is only visible in each of the test region(T), the test is invalid. Another test should be run to re-evaluate the specimen. Please contact the distributor or the store, where you bought the product, with the lot number.

#### Note:

1. There is no meaning attributed to line color intensity or width.
2. The abbreviations of drug names in the product images are just for demonstration. Please refer to the labels on the packaging and the prints on the test strip to identify which drugs are included in your test.

### QUALITY CONTROL

Users should follow the appropriate federal state, and local guidelines concerning the frequency of assaying external quality control materials. Though there is an internal procedural control line in the test device of control region, the use of external controls is strongly recommended as good laboratory testing practice to confirm the test procedure and to verify proper test performance. Positive and negative controls should give the expected results. When testing the positive and negative controls, the same assay procedure should be adopted.

### LIMITATIONS

1. This test has been developed for testing urine samples only. The performance of this test using other specimens has not been substantiated.
2. Adulterated urine samples may produce erroneous results. Strong oxidizing agents such as bleach (hypochlorite) can oxidize drug analyses. If a sample is suspected of being adulterated, obtain a new sample.
3. This test is a qualitative screening assay. It is not designed to determine the quantitative concentration of drugs or the level of intoxication.
4. It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
5. A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
6. The test result does not distinguish between drugs of abuse and certain medicines.
7. A positive result might be obtained from certain foods or food supplements.

### MEANING OF SYMBOLS ON PACKAGE

Keep away from sunlight	For in vitro diagnostic use only	Tests per kit
Store between 4°C-30°C (40°F-86°F)	Consult instructions for use	LOT Batch code
Keep dry	CE marking of conformity	Expiry date
Do not reuse	Authorized/representative in the European Community	REF Catalog number
Manufacturer	Recyclable	

### QUESTIONS?

If you have any question, please contact us via email [service@healthcare-manager.com](mailto:service@healthcare-manager.com). To learn more, please visit us at: [www.healthcare-manager.com](http://www.healthcare-manager.com). Distributed by Easy Healthcare Corporation, 360 Shore Dr., Burr Ridge, IL USA 60527

### QUESTIONS AND ANSWERS

**The Drug Line is lighter than the Control Line. Does this mean some drug is present?**  
No. Any line next to the word Drug or the drug abbreviation (depending on the test you have purchased), no matter how dark or light, is considered a Negative result and is not detected.

**How soon can I read my results?**  
Certain over-the-counter medications or prescription drugs may cross-react with the Easy@Home® Drug Test and cause a false positive result.

**What are any factors that can affect the test result?**  
Certain over-the-counter medications or prescription drugs may cross-react with the Easy@Home® Drug Test and cause a false positive result.

**How soon after taking a drug can you detect it in urine with a Easy@Home® Drug Test?**  
Most drugs can be detected in urine with a Easy@Home® Drug Test within a few hours after taking the drug; however this can vary depending on the type of drug taken, the amount taken, the frequency of use, and the metabolism of the person being tested.

**How long can a drug be detected in urine?**  
Each drug is cleared by the body at different rates. Some drugs, for example marijuana, can stay in the body for up to several weeks after use.

**What cut-off levels do Easy@Home® Drug Tests use for detecting drugs in urine?**  
The cut-off level for each drug varies (depending on the type of drug) and is measured in nanograms/ng/ml.

**Although Easy@Home® Drug Tests are designed to detect a very small amount of a drug in urine, if the amount is below the established cut-off level, you may test negative for that drug even though you may have taken the drug.**

**How soon after taking a drug can you detect it in urine with a Easy@Home® Drug Test?**  
Most drugs can be detected in urine with a Easy@Home® Drug Test within a few hours after taking the drug; however this can vary depending on the type of drug taken, the amount taken, the frequency of use, and the metabolism of the person being tested.

**How long can a drug be detected in urine?**  
Each drug is cleared by the body at different rates. Some drugs, for example marijuana, can stay in the body for up to several weeks after use.

**What are any factors that can affect the test result?**  
Certain over-the-counter medications or prescription drugs may cross-react with the Easy@Home® Drug Test and cause a false positive result.

**How soon after taking a drug can you detect it in urine with a Easy@Home® Drug Test?**  
Most drugs can be detected in urine with a Easy@Home® Drug Test within a few hours after taking the drug; however this can vary depending on the type of drug taken, the amount taken, the frequency of use, and the metabolism of the person being tested.

**How long can a drug be detected in urine?**  
Each drug is cleared by the body at different rates. Some drugs, for example marijuana, can stay in the body for up to several weeks after use.

**What are any factors that can affect the test result?**  
Certain over-the-counter medications or prescription drugs may cross-react with the Easy@Home® Drug Test and cause a false positive result.

**How soon after taking a drug can you detect it in urine with a Easy@Home® Drug Test?**  
Most drugs can be detected in urine with a Easy@Home® Drug Test within a few hours after taking the drug; however this can vary depending on the type of drug taken, the amount taken, the frequency of use, and the metabolism of the person being tested.

**How long can a drug be detected in urine?**  
Each drug is cleared by the body at different rates. Some drugs, for example marijuana, can stay in the body for up to several weeks after use.

**What are any factors that can affect the test result?**  
Certain over-the-counter medications or prescription drugs may cross-react with the Easy@Home® Drug Test and cause a false positive result.

**How soon after taking a drug can you detect it in urine with a Easy@Home® Drug Test?**  
Most drugs can be detected in urine with a Easy@Home® Drug Test within a few hours after taking the drug; however this can vary depending on the type of drug taken, the amount taken, the frequency of use, and the metabolism of the person being tested.

**How long can a drug be detected in urine?**  
Each drug is cleared by the body at different rates. Some drugs, for example marijuana, can stay in the body for up to several weeks after use.

**What are any factors that can affect the test result?**  
Certain over-the-counter medications or prescription drugs may cross-react with the Easy@Home® Drug Test and cause a false positive result.

**How soon after taking a drug can you detect it in urine with a Easy@Home® Drug Test?**  
Most drugs can be detected in urine with a Easy@Home® Drug Test within a few hours after taking the drug; however this can vary depending on the type of drug taken, the amount taken, the frequency of use, and the metabolism of the person being tested.

**How long can a drug be detected in urine?**  
Each drug is cleared by the body at different rates. Some drugs, for example marijuana, can stay in the body for up to several weeks after use.

**What are any factors that can affect the test result?**  
Certain over-the-counter medications or prescription drugs may cross-react with the Easy@Home® Drug Test and cause a false positive result.

**How soon after taking a drug can you detect it in urine with a Easy@Home® Drug Test?**  
Most drugs can be detected in urine with a Easy@Home® Drug Test within a few hours after taking the drug; however this can vary depending on the type of drug taken, the amount taken, the frequency of use, and the metabolism of the person being tested.

**How long can a drug be detected in urine?**  
Each drug is cleared by the body at different rates. Some drugs, for example marijuana, can stay in the body for up to several weeks after use.

**What are any factors that can affect the test result?**  
Certain over-the-counter medications or prescription drugs may cross-react with the Easy@Home® Drug Test and cause a false positive result.

**How soon after taking a drug can you detect it in urine with a Easy@Home® Drug Test?**  
Most drugs can be detected in urine with a Easy@Home® Drug Test within a few hours after taking the drug; however this can vary depending on the type of drug taken, the amount taken, the frequency of use, and the metabolism of the person being tested.

**How long can a drug be detected in urine?**  
Each drug is cleared by the body at different rates. Some drugs, for example marijuana, can stay in the body for up to several weeks after use.

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