
Smart Basal Thermometer BT-A31

User Manual

Guidance and manufacturer's declaration – electromagnetic emission		
The BT-A31 is intended for use in the electromagnetic environment specified below. The customer or the user of the BT-A31 should assure that it is used in such an environment.		
EMISSION TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
RF emissions CISPR 11	Group 1	The BT-A31 use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The BT-A31 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

CONTENTS	
1 INTENDED USE	01
2 CONTENTS	01
3 WHY BT-A31?	01
4 SET UP AND USE	01
5 DATA STORAGE	02
6 CHANGE THE BATTERY	03
7 TROUBLESHOOTING	03
8 TECHNICAL SPECIFICATIONS	04
9 MAINTENANCE & CLEANING	04
10 WARRANTY	05
11 WORKING & STORAGE CONDITIONS	05
12 PRECAUTIONS	06
13 EMC/FCC INFORMATION	06
14 WARNINGS	12
15 DISPOSAL OF THE DEVICE	12
16 EXPLANATION OF SYMBOLS	13


Guidance and manufacturer's declaration – electromagnetic immunity			
The BT-A31 is intended for use in the electromagnetic environment specified below. The customer or the user of BT-A31 should assure that it is used in such an environment.			
IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±8 kV air	±8 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: UT is the a.c. mains voltage prior to application of the test level.			

1. INTENDED USE	
The digital thermometer model No. BT-A31 is intended for the measurement and monitoring of human body temperature by doctor or consumers in the hospital or at home.	
2. CONTENTS	
1 Thermometer 1 User Manual	
3. WHY BT-A31?	
• Upon setting up, data auto sync to free app Premom. • Large and backlit LCD display helps you read results easily. • 1/100th High-precision with memory recall of last 30 readings.	
4. SET UP AND USE	
① Refer the quick guide of Premom App to set up the App before BBT measurement.	
② Take your temperature	

5. DATA STORAGE	
Memory stores up to 30 measurements on the device. While the thermometer is off double-click the On/Off button to review the stored data. You can exit the memory display mode by holding down the On/Off button for 3 seconds or it will shut off automatically after 30 seconds.	

6. CHANGE THE BATTERY	
Please replace the battery when the low-battery symbol appears on the right of the LCD display. Remove the battery cover and replace it with a battery of the same battery type CR2032H. Make sure the "+" sign faces up. Always dispose of used batteries according to local laws and regulations.	
① Symptoms: Unstable body temperature readings	
Possible Solutions: Ensure the thermometer probe is placed firmly underneath your tongue, close your lips, breathe normally and do not speak, eat, or drink.	
② Symptoms: Abnormal high BBT	
Possible Solutions: Take measurement immediately after waking up, limit physical activities before and during measurement.	
Note:	
1. If the problem persists, please contact customer support or consult your doctor if concerned.	
2. If you have any problems with this device, such as setting up, maintaining or using, please contact with Easy@Home customer service. Don't open or repair the device by yourself.	

8. TECHNICAL SPECIFICATIONS	
Product name	Electronic Clinical Thermometer
Model	BT-A31
Size	13.2cm x 2.92cm x 1.39cm
Transmission distance	9m
Display range	32.00°C - 42.99°C
Accuracy	±0.1°C(from 35.0°C to 39.0°C) ±0.2°C(from 32.0°C to 35.0°C, and from 39.0°C to 42.99°C)
Battery	CR2032H
Receiver	Apple Phone: iPhone5+ IOS8.0+ Android Phone: Android 4.4+
Battery life	about 3 months
Operating conditions	Temperature: 5°C to 40°C
Weight	About 20g (without battery)

9. MAINTENANCE & CLEANING	
1. Clean the tip of the thermometer with a solution of 75% rubbing alcohol and 25% water before and after every use.	
2. Do not use cleaning agents other than alcohol or water to clean the thermometer as it may damage or decrease the lifetime of the product and / or present safety risks.	
3. While you can clean the sensor tip of the thermometer with water, it's important to avoid fully submerging the device in water or any other liquid. Doing so can result in damage to the thermometer, potentially leading to inaccurate readings and posing a safety risk.	

14. WARNINGS	
1. No servicing/maintenance while the thermometer is in use.	
2. Not for use in an OXYGEN RICH ENVIRONMENT.	
3. Before every use, check the device. Do not use the device or an electrode if it is damaged in any way. The continuous use of a damaged unit may cause injury, improper results, or serious danger.	
4. Don't open or repair the device by yourself.	
5. The housing materials of the thermometer are ABS and TPE. Be careful of potential allergic reactions to these materials.	
6. The typical service life of the new and unused batteries is 100 hours for continuous operation.	
7. Protection against electric shock: Internally powered ME equipment.	
8. Protection against harmful ingress of water or particulate matter: IP22. However, for the sake of accuracy and safety, please refrain from fully immersing the device in water or any other liquid.	
9. Keep new and used batteries out of the reach of children, and seek immediate medical attention if a battery is ingested, and to follow any other consensus medical advice.	
10. Choking Hazard-Small parts not for children under 3 years or any individuals who have a tendency to place inedible object in their mouths.	

16. EXPLANATION OF SYMBOLS	
	CE marking of conformity
	Type BF applied part
	Manufacturer
	Authorized representative in the European Community
	WEEE (Waste Electrical and Electronic Equipment)
	Manufacturing date
	Batch code
	Keep away from sunlight
	Ingress protection rating
	Caution
	Humidity limitation of 0% ~ 85%
	Temperature limit of 20°C-55°C
	Recyclable
	Keep out of reach of children


www.premom.com
Questions or comments?
E-mail: support@premom.com
Distributed by: Easy Healthcare Corporation
360 Shore Dr., Burr Ridge, IL USA 60527
Made in China

12. PRECAUTIONS	
1. Do not bite, bend, drop or take apart the thermometer.	
2. Keep the device out of the reach of children/pets to avoid inhalation or swallowing of small parts.	
3. Do not expose it to direct sunlight, high temperature and moisture.	
4. Not intended to be sterilized. Prevent saliva or cleaning solution from penetrating the display window.	
5. Clean the thermometer after or before use with soft cloth and cleaning solution listed below and store in the plastic case. Recommend to wipe with wet water soaked cloth for home use.	
6. If LCD shows the symbol or LCD shows unclear, it means the battery has run out. Please replace the battery within 10 minutes.	

13. EMC/FCC INFORMATION	
1. This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile RF communications equipment.	
2. Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.	
3. Caution: This unit has been thoroughly tested and inspected to assure proper performance and operation!	
4. Caution: This machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.	

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