
### BianStone & Graphene Far Infrared Heating Therapy MagicPad

#### Product Introduction


This revolutionary MagicPad adopts a brand-new heat energy injection technology, and it only takes 15 minutes to bring you an unparalleled feeling of total soothing and relaxation in an instant. Human body safety 24V/3A low-voltage input can heat up extremely quickly, emit far infrared rays and negative ions into the body, generate the same frequency resonance with cells, activate cells, promote blood circulation, relieve pain, improve deficiency-cold constitution and sleep quality, warm compress to reduce inflammation, restore physical fitness, stabilize emotions, and help the body absorb nutrients. Long-term use can promote metabolism and enhance human immunity. Using high-purity Graphene heating cloth is used to heat, not by wire, the harmful electromagnetic wave is close to 0, waterproof, fireproof, and electric shock proof. It can be used as a cushion or foot pad, as well as hot compresses for shoulders, stomachs, and waist. Whether at home or at work, it is your good friends/enjoy pure natural BianStone (scraping stone) & Graphene far-infrared warming health care.

- #### Eight characteristics
- Natural premium BianStone
  - Graphene heating Not by wire
  - Far-infrared wavelength 6-15 microns
  - Anti-burning Anti-leakage
  - Negative ions
  - Multifunctional belt
  - Sustainable temperature 20-60°C
  - 3D shockproof

1

### Product usage diagram


#### Product connection diagram:

- Connect the four-pin connector (male) with the thermostat's four-pin connector (female) according to the limit slot and arrow to confirm the built connection.
- Connect the DC connector (female) at the other end of the thermostat and the DC connector (male) of the power supply end to the limit slot attached to it and rotate and fix it;
- Insert the end of the matching power cord into the power connector.


2

### ④ Multifunctional belt

3

### Controller use and troubleshooting

- After the power is turned on, long press the key to turn on or off. The default temperature "30°C" is set to "8" hours when it is used for the first time.
- The icon is for temperature adjustment. Press or to adjust for temperature rise and fall. Note: After the temperature is set, the displayed temperature will automatically jump to the actual sensing temperature.
- The icon is for timing adjustment. Press or for timing adjustment. Note: The timing duration is fixedly set to 2 hours, 4 hours, 6 hours, 8 hours, and 12 hours, and it will automatically stop shutting down after reaching the setting.
- After the temperature is set, the icon flashes to indicate heating is in progress, and when the temperature is reached, the icon disappears to stop heating.
- The controller has a memory function, after the last temperature and time adjustment, the default last adjustment parameters will be used again;

mode	set up	Automatic heating time
A.Warm therapy 1	20-45°C	Can be set manually for 2-12 hours
B.Warm therapy 2	46-50°C	Heating for 4 hours, automatically return to 37°C
C.Hot therapy	51-60°C	Heating for 4 hours, automatically return to 37°C

#### 6. Thermostat display troubleshooting

display	cause of issue	Troubleshooting
E1	1. The thermostat is not connected to the product or the connection is in poor contact. 2. The temperature sensor probe fails.	1. Check the connection line, whether the connection is normal, and it needs to be re-wired. 2. If the product is still abnormal, please contact after-sales service.
EH	High temperature error or temperature out of protection	Turn off the power, turn on the product again after cooling to confirm, if it shows that the product is abnormal, please contact after-sales

4

### The product soothes and improves the following people

- Insomnia, emotional imbalance, long-term fatigue, deficiency and cold constitution;
- Women with physical pains, cold hands and feet, edema and obesity, sub-health;
- Back pain, shoulder and neck pain, numbness of hands and feet;
- Sports injuries.

#### Physiotherapy mode

- Daily Warm Compress Therapy:** Set the thermostat to 30~40°C, for daily use. It is suitable for those with deficient constitution, poor body microcirculation, long-term fatigue, rheumatism pain, cold hands and feet, weakened immunity. It is recommended to use heating for 6 to 10 hours/time.
- Warm Conditioning Therapy:** The thermostat is set to 41~49°C. For daily body conditioning. It is suitable for resisting chronic diseases, long-term fatigue, rheumatic pain, dysmenorrhea, uterine cold, blood stasis, cold hands and feet, poor body microcirculation, decreased immunity, edema, obesity, sub-health; hemorrhoids, male prostatitis and other symptoms to restore health. It is recommended to use 2 to 4 hours / time for heating.
- Hot Therapy:** The thermostat is set to 50°C~60°C, it can be used for heat therapy when necessary. It is used for sports injuries, backache and bone pain, shoulder and neck pain, and numbness in hands and feet. It can accelerate blood circulation, perspiration and detoxification, and effectively relieve pain. Recommended use of heating: 1 hour / time.
- Natural Therapy:** The thermostat power heating switch is turned off (room temperature), natural care on hot days or normal indoor temperature. It feels cool and comfortable when used in Summer; in normal long-term use, the body can absorb the far infrared rays that are continuously and naturally emitted by the BianStone, which is beneficial to the body, which is beneficial to maintain and improve the normal operation of the body's microcirculation system. Recommendations for use: long-term.

#### Specifications

Model : SP450P  
Product name : BianStone & Graphene Far Infrared Heating Therapy MagicPad  
Size : 45cm x 45cm Weight : 1.1 KG  
BianStones quantity : 59 PCS  
Therapy mode : Room temperature natural therapy, warm therapy, hot therapy  
Adapter input voltage : 100-240V-50/60hz 2.0A  
Adapter output voltage/current : 24V 2.5A  
MagicPad power : 72W max  
MagicPad Current and Voltage : Human body safety low voltage max 24V/3A  
Heating device : Graphene heating film  
Adjustable range of temperature controller : 20-60°C  
The temperature controller can set the timing : 2-12 hours  
Heating temperature : 20-60°C up and down button adjustment  
Accessories : 3D shock net, Density mesh cloth, Bamboo charcoal fiber cotton (freeproof)

5

### Package Contents

- \* MagicPad x1
- \* Manual book x1
- \* Temperature controller x1
- \* Power adapter x1

#### Warranty description

Zikko provides quality warranty services for hardware products and accessories in the packaging and defects in materials and workmanship. The warranty period of the Dr.Rock MagicPad (including temperature controller) is one year from the date of purchase of the product, and the included storage bag is guaranteed for 6 months. Damages caused by normal wear or aging, accidents, improper use, alterations, etc. are not covered by this warranty. How to obtain after-sales service: Please call the Zikko authorized brand or go to the authorized retailer's purchase place or online customer service for after-sales service, and please provide the relevant purchase credentials (invoice or purchase record).

#### Product maintenance methods

You can use a wet towel to wipe the contaminated part when it is not powered on, or use a disinfectant wet towel to disinfect it regularly, and then turn it on for heating after drying.

#### ⚠ Precautions

- Read the manual carefully before use and use it in accordance with the operating instructions.
- Don't let it be used alone by people who cannot take care of themselves, infants or children who are not sensitive to heat.
- It is forbidden to use in water and humid environment.
- Do not dry clean or wash. You can wipe the contaminated part with a wet towel when the power is not turned on, and then turn on the power to heat it after drying.
- Do not insert sharp objects into this product.
- Do not place bedding and other items in a small area to cover the heating area when the product is heated to prevent local overheating and causing malfunctions.
- It is forbidden to place flammable materials such as fireworks and lighters on the surface of the product to prevent fire.
- Please confirm whether the power supply is cut off before going out.
- Using the heating function, the temperature can be adjusted according to the tolerance of the skin. In order to prevent skin burns. If you feel unwell, stop using it in time.

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