breo.iBack2

Massage Pillow
Coussin de Massage
Massagekissen

breo

Product Configurations

Inside the Package

massager

adapter

user manual

Product features

Warm Compress

Light & Portable

Multi-function

Kneading Massage

Preparations

Please make sure that the device and accessories are complete and the power adapter is not damaged before use.

Operations

Power on

Press and hold the power button for 2 seconds to turn on the massager.

Heating function

The heating function is deactivated by default, and can be turned on or off by short pressing the power button.

Automatic shut-off

The massager automatically turns off after a 15 minutes' session.

Power off

Press and hold the power button for 2 seconds to turn off the massager.

Note

The massager functions only when it is connected to the mains power. It is not equipped with rechargeable batteries, avoid trying to charge it.

Precautions

Please read this User Manual carefully before using this device to ensure safe and proper use.

The company shall not be responsible for any consequences caused by improper use or non-compliance with the relevant instructions of this device.

1. Do not use the device to massage any area that has been operated on before. Please follow your doctor's advice if you need to use the device.
2. Do not use the device to massage any area that is suffering from any sort of discomfort such as skin injury, swelling, inflammation and bruising.
3. People whose movements are slow, or who are physically or mentally disabled, or mentally impaired should not use it without authorization. Please follow your doctor's advice if you need to use the device.
4. This device has been designed as an electrical healthcare product and is not intended for diagnostic, therapeutic or other commercial purposes.
5. Children under 8 years of age are forbidden to use the product; children between 8-12 years of age must use the product under the guidance and care of their parents or other persons with supervision capability.
6. Do not use the device while driving. Please abide by the safe driving rules.
7. Please only use the power adapter provided by our company and ensure electrical safety.
8. To clean the device exterior, please remove the cloth cover and wash with water.
9. Do not hit, throw, or use sharp objects to scratch or poke the product.
10. Minors must only use the product under guidance from their parents or persons with supervision authority.
11. If there are any abnormalities during use, please cut off power immediately and contact our company. Do not disassemble and repair it on your own.

12. Please stop using the product if you experience any physical discomfort or abnormalities during use.

13. The device has a heating surface. Those who are not sensitive to heat must be cautious when using the device.

Troubleshooting

If a fault occurs during the operation of the device, please check the following precaution.

If the error cannot be fixed or the device still does not work properly, please cut off the power and contact Breo to repair the device. Do not disassemble or try to repair the device yourself.

If you have any issues, questions or feedback, please send an email to support@breocare.com or call us at 866-546-3587. We are always happy to help.

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Dimensions: Unité principale: L370 mm × l220 mm × H110 mm

Fabriqué en Chine

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Shenzhen Breo Technology Co., Ltd
19/F, Hisense Southern Building,
No.1777 Chuangye Road, Yuehai Street, Nanshan District,
Shenzhen 518054, China