

Product Information

Product Description

This product is a multi-polar RF home beauty device ("RF device") designed to improve facial problems experienced by women using radio frequency ("RF"), electrical muscle stimulation ("EMS"), and triple effect phototherapy.

Product Specifications

- Product name: AMIRO R1 Lift Facial RF Skin Tightening Device
- Product model: ARF001
- Product size: 160×38×38mm
- Size of electrode contact surface: 32×32mm
- Rated input: 5V = 2A
- Battery capacity: 2600mAh
- Country of origin: China

Product List

AMIRO R1 Lift Facial RF Skin Tightening Device x 1

AMIRO BEAUTY Moisturizing Gel x 1


Charging cable x1

User manual x1

Instruction

Product illustration

1. Electrodes (6)
2. Temperature sensor
3. LED (internal)
4. Mode indicator
 - ☀ Eye care mode
 - ☾ Full face anti-aging mode
5. Power/mode button
6. Intensity "+" button
7. Intensity "-" button
8. Intensity/power indicator
9. Type-C charging port



Operating procedure

1. Charging

- Insert the Type-C connector of the charging cable into the device's charging port.
 - The device is fully charged when all 5 intensity/power indicators are lit.
- **The device cannot be used when charging.
- **Lay the device down when charging to prevent it tipping over.

2. Powering on and off

- When powered on for the first time, the device will be in default eye care mode. The device will be in the last-used mode when it is powered on again.
- The device will automatically turn off if no operation and no contact with skin are made within 120 seconds after being powered on.
- Long press the power/mode button for 2s to power off the device.

3. Switching modes

- When power is on, short press the power/mode button to switch between [eye care mode] indicated by "▲" and [full face anti-aging mode] indicated by "▼".

4. Selecting intensity

- The default intensity is 1 (1 light on) when first powered on; it remembers the last intensity used when it is next powered on.
- After powering on the device, short press the "+" button to increase intensity and the "-" button to decrease intensity.
- The highest intensity is 5 (5 lights on), and the lowest is 1 (1 light on).

5. Power indicator

- The device does not display remaining power when powered on or used.
- Press the [power/mode button] to check the [intensity/power] indicator for remaining power when powered off. 5 indicators indicate the highest power level and 1 indicator indicates the lowest.

6. Low battery and power shortage warnings

- In powered-on or standby states, the corresponding indicator will keep flashing when power is low. Please charge regularly to ensure smooth use.
- When powered off, the power is low if the device cannot be turned on by pressing the power/mode button and the Level 1 indicator flashes 5 times before going out. Please charge immediately to ensure normal operation.

Care advice

- **[Eye care mode]: It is recommended to use 2-3 times per week, in combination with AMIRO BEAUTY Gel or eye mask. Usage sequence: Left eye area and right eye area for 2 minutes each, 4 minutes in total.**
- **[Full face anti-aging mode]: It is recommended to use 2-3 times per week, in combination with AMIRO BEAUTY Gel or facial and neck mask. Usage sequence: V-shaped face mask for 2 minutes, left cheek for 1 minute, right cheek for 1 minute, forehead for 1 minute, 5 minutes.**

Care procedure

Step 1: Makeup removal and cleansing

- Remove makeup and cleanse your face to ensure your skin is clean.

Step 2: Applying the gel

- Apply some recommended gel to the center of the RF electrode or the face (When you feel the resistance of the device on the face

or skin dryness, please replenish the gel promptly)

[Eye care mode] for 4 minutes. Start the RF device and select [eye care mode] “▲”

- Apply AMIRO BEAUTY Gel or eye mask around the eyes
- Place the RF device on the left eye area and glide. When using the device for the first time or after not using it for a long period, start with the lowest intensity to build tolerance (see illustrations below)
- After the second brief vibration, move the device onto the right cheek



[Full face anti-aging mode] The whole care process is about 5 minutes. Turn on the radio frequency instrument, and select [full face anti-aging mode] "▼."

- Apply AMIRO BEAUTY Gel evenly on the entire face or place a face and neck mask on the jawline area
- Place the RF device on the V-shaped face mask area and glide. When using the device for the first time or after not using it for a long period, start with the lowest intensity to build tolerance (see illustrations below)
- Move the device to the left cheek area after the first brief vibration, to the right cheek after the second brief vibration and to the forehead after the third short vibration.
- When the device vibrates for 1 second and is powered off, the treatment mode is finished. Switch it on again if you need to continue using it.



Step 4: Cleansing and hydrating

- After the care routine, wash your face with clean water and apply skin care products if needed.

Step 5: Device maintenance

- After treatment, wipe off any gel remaining on the RF device with a slightly damp soft cloth, then dry it with a dry cloth before putting it back in a box or storage bag.

Treatment techniques (illustrations and introduction)



- **[Eye care mode]:** 4 minutes, 2-3 times per week

- **[Full face anti-aging mode]:** 5 minutes, 2-3 times per week

Tips:

- Place the device against the face. The care begins after a brief vibration. Move to a different area after each vibration.
- When the device vibrates for 1 second and is powered off, the treatment mode is finished.
- Restart and short-press the mode button to switch to another mode and repeat the above steps. Use no more than 30 minutes per day.
- During use, please move the RF device slowly for the best effect and do not stop.

Safety precautions

Tips:

Carefully read all [Safety Precautions] before using the product and adhere to all guidelines in this handbook to prevent adverse effects.

[Use area]:



Operation precautions

- The recommended ambient temperature for this product is 10°C-30°C. Ensure the temperature in the charging environment is between 10-30°C.
- When using the device for the first time or if a long time has passed since the last use, first charge the device and make sure all buttons and indicators are normal.
- Set intensity to Level 1 when using the device for the first time or before your skin gets used to this device
- If you experience abnormal skin conditions, stop using the device immediately and contact our after-sales service or consult a professional doctor.
- When using the device, ensure the RF electrodes are pressed against the skin. Do not glide the device quickly over your skin. The recommended speed of movement is 2cm/s; if it is too hot, you may increase the speed slightly.
- Do not allow the device to remain in the same position over your skin for long to prevent burns. Glide the device over your skin.
- Do not look directly at the LED for a prolonged period to prevent damage to eyes.
- Do not apply the RF electrodes to the eyelids or around the eyes to prevent damage to eyes
- Differences in individual skin sensitivity may lead to differences in temperature perception.
- When the RF device is connected to the power supply for charging, it cannot be turned on or work and all the keys are invalid.
- Remove the battery before discarding the device, and dispose of the battery safely. Battery removal steps: Ensure that the product is powered off, and then you may use tools to assist you in detaching the bottom cover and removing the battery inside. The battery should be recycled as required.

Precautions

- This product is for people over 18 years old; please do not use it on minors; please keep it in a place out of reach of children.
- Do not use the device on pets as it may cause injury to them.
- Do not use the product for purposes other than those specified in the manual to avoid personal injury or mechanical failure.
- People with physical or mental disabilities should not use this device as it may result in accident or bodily injury.
- Do not use the device if it is malfunctioning as it may result in bodily injury or danger.
- Do not use nearly flammable substances (such as alcohol, thinner, gasoline, nail polish remover or styling spray), which may result in fire, electric shock or bodily injury.
- Do not use the device together with flammable cosmetics.
- Do not charge, use or store the product near a fire, heat source, or high temperature area, as it may cause the device to overheat, catch fire, or burst.
- Do not use near water sources, in bathrooms or while bathing, as it may result in electric shock, short circuit, or fire.
- Please avoid dropping the device and the base or collisions to prevent malfunctions or injuries due to the fall. Do not use the device if it has been damaged by a fall or collision.
- Avoid contact of the electrodes with metals as it may cause device malfunction and personal injury.
- Do not use corrosive solvents (alcohol, thinner, gasoline, nail polish remover, etc.) to wipe the device, as it may damage its surface
- Turn the power off when the device is not in use or being cleaned to prevent electric shock, short circuit, or fire.
- Please do not disassemble or modify the device, as this may cause fire, electric shock, or physical injury; if the device needs repair, please contact our after-sales service.
- If any abnormalities occur during charging, please unplug the charging cable immediately. Please cease operating the equipment immediately if there is smoke, noise, or an odor.

The device is unsuitable for the following groups of people:

People with any of the following conditions may not use this product as it could result in accidents or skin or bodily injury:

- Menstruation, pregnancy or breastfeeding
- Puberty or adolescence
- Metal implants around the application areas

- Diseases such as cancer (especially skin cancer and pre-cancerous moles), acute diseases, tuberculosis, febrile diseases, infectious diseases, blood diseases, visceral diseases, skin diseases, allergic dermatitis, etc.
- For those with allergies, skin conditions, allergic dermatitis, sensitive skin, etc., using this product may cause discomfort or adverse effects. In case of adverse effects, immediately cease its use
- Impaired immune system due to immunosuppressive diseases such as AIDS or taking immunosuppressive drugs
- Consult a doctor if you have had permanent plastic fillers within the application area, or you have temporary dermal fillers or Botox implanted within the application area
- Any disease within the application area due to thermal irritation, such as recurrent herpes simplex virus
- Possible blemishes within the application area cause by physical irritation
- Allergic, sensitive or susceptible skin
- Fatigue, illness, under the influence of alcohol or sleeping pills
- Dilated capillaries caused by prolonged use of steroids, hormones, and other drugs, or liver dysfunction
- People using these medical devices: Assistive devices implanted in the body such as artificial lungs, pacemakers, internal defibrillators, or portable devices such as electrocardiographs.

Where the device should not be applied:

Do not use this product on the following body parts or if the following conditions exist:

- Throat nodules, eyeballs, eyelids, oral cavity, private parts, mucous membranes, etc.
- Severe acne in the application area
- Purulent and inflamed skin, previous scarring or abnormal wound healing
- If you have undergone surgeries such as laser peeling or chemical peeling in the last 3 months or have healing wounds within the application area
- Skin inflammation or redness due to sunburn
- Itchiness or pathological pigmentation on the skin
- Protruding moles or viral warts on the skin
- Dermatological conditions such as rashes, inflammation, itchiness, etc.
- Areas with damaged skin integrity (e.g., burns, blisters, cuts, wounds, active skin disease or infection).

Note: The above list of safety precautions is not exhaustive. If you are unsure whether you can use this product, please consult a doctor or contact us.

Warranty

- AMIRO provides 1 year limited warranty.
- For more information about the warranty policy, please refer to AMIRO's official website www.amirobeauty.com.
- Except as prohibited by laws or otherwise promised by AMIRO, the after-sales services shall be limited to the country or region of the original purchase. Under the consumer warranty, to the fullest extent permitted by law, AMIRO, will, at its discretion, repair, replace or refund your product.
- Normal wear and tear, force majeure, abuse or damage caused by the user's negligence or fault are not warranted. If in doubt please contact the relevant person as AMIRO may identify. Products which were not duly imported and/or were not duly manufactured by AMIRO and/or were not duly acquired from AMIRO or a AMIRO's official seller are not covered by the present warranties. As per applicable law you may benefit from warranties from the non-official retailer who sold the product. Therefore, AMIRO invites you to contact the retailer from whom you purchased the product.

Malfunction and maintenance

Table of common malfunctions

Problem	Possible cause	Solution
Start failure	Dead battery	Charge fully before use
	End of battery life	Contact the manufacturer for repair or replacement
Automatic power off	No operation for more than 120 seconds	Restart the device
	Dead battery	Charge fully before use
	RF device electrodes not tightly pressed against the skin	Keep RF device electrodes tightly pressed against the skin
	Inappropriate gel	AMIRO BEAUTY gel is recommended
	Insufficient gel	Replenish gel on the face
Burning sensations of the face	High intensity	Lower intensity
	RF device stays in one area	Slightly speed up RF device movement
	Insufficient gel	Replenish gel on the face
Low heat	Low intensity	Increase intensity
	Fast RF movement	Slightly reduce RF device movement
	Inappropriate mode	Around the eyes special care mode RF energy for the skin around the eyes designed for the whole face lifting mode to the current sense
Tingling caused by current	Device not pressed against the skin	Press the device against the skin
	High mode intensity	Lower intensity
Redness of skin after use	Caused by RF heat	It normally recovers soon afterwards; please seek medical attention if you don't recover or have other discomfort.



Warning:
Non-ionizing radiation



Read user manual