

Q1:How to use this thermometer?

Step 1: Press the power button to turn it on.

Step 2: Wait for all the display symbols to appear.

Step 3: The thermometer will display the last recorded measurement from memory.

Step 4: Allow approximately 10-15 seconds for the measurement to complete. The thermometer will indicate completion with a series of "Bi-Bi-Bi-Bi" beeps.

Q2: Are rectal thermometers better for babies than other thermometers?

Yes, rectal thermometers are preferred for babies. Other options like armpit and oral thermometers can be uncomfortable and difficult to use for infants and some older adults. Pediatricians recommend rectal thermometers for children under 3 months as they are more reliable and accurate. They measure internal body temperature, reducing the influence of external factors and providing more precise readings.

Q3: How long does it take to obtain a temperature reading with this thermometer?

The test time for this thermometer is relatively quick, typically ranging from 10 to 15 seconds. This allows for a swift and efficient measurement process.

Q4: Does the thermometer's probe have the capability to be bent?

Yes, the flexible probe allows for effortless, secure, and gentle measurement.

We hope these tips were useful. If you have any inquiries about our products or services feel free to contact us via email, call our customer service at 855.822.6999 M-F 9am-5pm CT or visit us at [healthcare-manager.com](https://www.healthcare-manager.com)