

Q1: How does the Dual Color Backlight and Fever Alert feature work?

The thermometer's Dual Color Backlight and Fever Alert feature indicates different temperature ranges. A normal temperature range is displayed with a green backlight and a 10-second-long consistent beep sound. However, a high temperature range is indicated by a red backlight and rapid beep sounds for 10 seconds.

Q2: Can I access previous temperature readings with this thermometer?

Yes, the thermometer has a Memory Display function. It can recall the last 30 temperature measurements, allowing you to view your previous readings.

Q3: How does the thermometer exit memory recall mode?

The thermometer will exit memory recall mode in two situations. Firstly, it will exit when it reaches the oldest recorded measurement in memory. Secondly, if you have not pressed any button for approximately 1 minute, it will automatically exit memory recall mode.

Q4: Does the digital oral thermometer come with an additional battery?

Yes, the packaging of the digital oral thermometer includes an additional battery. This provides you with a spare battery to replace the existing one when needed.

Q5: How do I replace the battery when the low battery icon appears?

When the low battery icon appears in the lower right corner of the display, it indicates that the battery needs to be replaced. To replace the battery, simply slide open the battery cover, insert the new battery with the positive (+) side facing up, and then replace the cover securely.

We hope these tips were useful. If you have any inquiries about our products or services feel free to contact us via email, call our customer service at 855.822.6999 M-F 9am-5pm CT or visit us at [healthcare-manager.com](https://www.healthcare-manager.com)