




Autoclavable Indicator Tape

Used to indicate that instruments or items have been thoroughly sterilised. To ensure proper steam sterilisation, the steam autoclave process must reach and maintain 121°C for 15 minutes, and if it's a higher temperature of 134°C, 3.5 minutes to 4 minutes is the standard time. Our TST strips or steam indicator strips will change colour from blue to black once the temperature has been reached. Available in 3 sizes.

CODE	SIZE	BUY 1
990670	13mm x 55metres	\$3.85/roll
990671	19mm x 55metres	\$4.85/roll
990672	25mm x 55metres	\$6.25/roll



TST Strip Class 4

BUY 1 @ \$15.95/box


200pcs/box

New waterproof design. Used to indicate that instruments or items have been thoroughly sterilised. To ensure proper steam sterilisation, the steam autoclave process must reach and maintain 134°C for 2 minutes. Our TST strips will change colour from green to black once the temperature, steam pressure and time has been reached.

Size: 100mm (L) x 18mm (W)

**FREE
SAMPLES!**

CODE
990649



**Each strip can be cut
to make two pieces.**

TST Strips Class 4

BUY 1 @ \$3.25/box

100pcs/box

Used to indicate that instruments or items have been thoroughly sterilised. To ensure proper steam sterilisation, the steam autoclave process must reach and maintain 134°C for 2 minutes. Our TST strips will change colour from blue to black once the temperature, steam pressure and time has been reached.

Size: 90mm (L) x 20mm (W)

**FREE
SAMPLES!**

CODE
990650



TST Indicator Strips Class 6, 8 min

BUY 1 @ \$39.95/box

200pcs/box

Used to indicate that instruments or items have been thoroughly sterilised. To ensure proper steam sterilisation, the steam autoclave process must reach and maintain 135°C for 8 minutes. Our TST indicator strips will change colour from purple to green once the temperature, steam pressure and time has been reached.

Size: 62mm (L) x 25mm (W)

**FREE
SAMPLES!**

CODE

990652



TST Indicator Strips Class 6, 3.5 min

BUY 1 @ \$39.95/box

200pcs/box

Used to indicate that instruments or items have been thoroughly sterilised. To ensure proper steam sterilisation, the steam autoclave process must reach and maintain 135°C for 3.5 minutes. Our TST indicator strips will change colour from purple to green once the temperature, steam pressure and time has been reached.

Size: 62mm (L) x 25mm (W)

**FREE
SAMPLES!**

CODE

990653



Green Card T, Bowie Dick Test

BUY 1 @ \$39.95/bag

15pcs/bag


Used to detect the presence of residual air in small pre-vacuum steam sterilisers operating at 135°C for 3.5 minutes. The card changes colour from purple to green when exposed to a steam sterilisation cycle with acceptable vacuum cycles.

Size: 85mm (L) x 85mm (W)

**FREE
SAMPLES!**

CODE

990643




Bowie Dick Test, Tray

BUY 1 @ \$32.95 each

A tray for the Green Card T, to be placed into the steam steriliser. Autoclavable up to 135°C.

CODE

990644




Self-sealing Sterilisation Pouches

200pcs/box

With internal steam indicator. Triple adhesive seal to protect contents. Blue tinted plastic window helps to detect puncture holes that may be caused by sharp or jagged instruments. Disposable.

**FREE
SAMPLES!**

CODE	INNER POUCH SIZE	BUY 1
990621	57mm x 130mm (2 1/4" x 4")	\$3.95/box
990614	50mm x 230mm (2" x 9")	\$5.55/box
990619	70mm x 230mm (2 3/4" x 9")	\$6.95/box
990620	90mm x 145mm (3 1/2" x 5 3/4")	\$6.95/box
990616	90mm x 230mm (3 1/2" x 9")	\$7.65/box
990617	135mm x 260mm (5 1/4" x 10")	\$14.25/box
990618	190mm x 330mm (7 1/2" x 13")	\$22.95/box



Sterilisation Pouch Dispenser-Medium

BUY 1 @ \$79.50 each

Ideal for storing sterilisation pouches of popular sizes. Pre-drilled holes for wall-mounting. Location A can store width smaller than 60mm. Location B can store width smaller than 60mm or TST Strips. Location C can store width smaller than 90mm. Location D can store width smaller than 90mm. Location E can store width smaller than 135mm. Size: 440mm(L) x 81mm(W) x 307mm(H). Note: Pouches not included. Cleaning Instructions: Use only soap and water on acrylic surface. Alcohol or ammonia-based disinfectants may cause cracks to appear.

CODE	COLOUR
990639	Clear



Sterilisation Roll

BUY 1 @ \$18.95/roll

Various sized reels allows for customized sterilization of a wide variety of items.

Size: 50mm x 200metres

CODE
990632

100%
RECYCLED PAPER



CSR Sterilisation Wrap

Conforms easily to various instruments and items of various shapes and sizes. Designed for use with ETO, steam and gamma irradiation sterilisation. These wraps provide an effective bacterial barrier to both airborne and waterborne contamination.

CODE	SIZE	QTY	BUY 1
990655	300mm x 300mm	100pcs/bag	\$6.85/bag
990656	380mm x 380mm	100pcs/bag	\$11.95/bag
990657	500mm x 500mm	100pcs/bag	\$21.65/bag
990658	610mm x 610mm	100pcs/bag	\$29.65/bag
990659	760mm x 760mm	50pcs/bag	\$24.25/bag



Sterile Surgical Drape

ETO Gas Sterilised. Made of polypropylene with lamination. Ideal for dental implant surgeries to reduce bacterial transfer and contamination of the surgical site. These drapes are highly absorbent. Clean-up time is greatly reduced.

CODE	SIZE	QTY	BUY 1
990720	450mm x 750mm	20pcs/box	\$13.95/box
990721	600mm x 600mm	20pcs/box	\$13.95/box
990722	900mm x 900mm	20pcs/box	\$19.95/box