VELPEAU[®]









[Product Name] : Knee Support

[Product Code]

[VP1206]

□ VP1206-S □ VP1206-M □ VP1206-L □ VP1206-XL □ VP1206-XXL [VP1206C] :

□ VP1206C-S □ VP1206C-M □ VP1206C-L □ VP1206C-XL □ VP1206C-XXL

[Product Specifications]:

	Size	Product code	Color	Thigh circumference (inch)	Left / Right
	S	VP1206-S	Grey	13"-14 1/2"	Exchangeable
		VP1206C-S	Black		
	М	VP1206-M	Grey	14 1/2"-16"	
		VP1206C-M	Black		
	L	VP1206-L	Grey	16"-17 3/4"	
		VP1206C-L	Black		
	XL	VP1206-XL	Grey	17 3/4"-20"	
		VP1206C-XL	Black		
	XXL	VP1206-XXL	Grey	20"-22 1/2"	
		VP1206C-XXL	Black		

Measure the circumference of thigh. Please measure 5 inches above the kneecap. Choose the appropriate product size according to the





www.velpeau.com

VELPEAU[®]



Shenzhen Weprotex Technology Product Co.,Ltd





www.velpeau.com -

[Product Features]:

- •The use of high-quality spandex, nylon, and other materials ensures excellent elasticity, providing a comfortable and secure fit that prevents slipping while wearing.
- ${\color{red} \bullet } \textbf{The knitted design follows ergonomic principles, applying balanced}$ pressure for effective adjunct therapy.
- •The seamless high-temperature compression technology ensures a comfortable wearing experience without friction on the skin..
- •The silicone pad fits tightly around the patella, providing stability, massage, shock absorption, and protection
- $\bullet \textbf{Enhances knee joint strength, reduces fatigue, improves endurance, and} \\$ offers optimal knee protection.
- •The soft and elastic metal spring support bars on both sides enhance knee joint explosive power during exercise.
- ·Features a stylish two-color design.
- ·Manufactured using advanced German machines, guaranteeing top-quality construction with no need for trimming.
- •Double silicone gel strips ensure secure placement during intense activities, preventing the product from slipping down.

[Product Compounds]: Nylon 45%, Silicone 20%, Spandex 18%, Polyester 9%, Stainless Steel 8%.

[Application]:

This product is designed to provide support and stabilization to injured knee areas, restricting movement and relieving pain. It serves as an adjunct therapy for partial damage, aiding in the recovery process.

[Precautions / Warnings] :

- •Please consult a professional for guidance before using this product.
- •Please use with caution if you have skin ulceration or open wounds. It is important to protect the affected area and avoid direct contact between the brace and the wound to prevent any potential discomfort or complications.

— www.velpeau.com —

【Transportation And Storage Conditions】:

The product can be transported using general means of transportation and stored at room temperature. However, it must be fireproof, moisture-proof, dustproof, rodent-proof, and insect-proof.

【 Washing 】:



W

Hand wash





X







 \boxtimes

Do not dry clean

[Term of Validity]: 5 years

[Production Date]: See the package

[Manufacturer]: Shenzhen Weprotex Technology Product Co., Ltd.

[Production Address]: 1710, Building 11, Phase II of Tian'an Yungu Industrial Park, Gangtou Community, Bantian Street, Longgang District, Shenzhen, China

[After-sale Service]: Shenzhen Weprotex Technology Product Co., Ltd.

[After-sale Service Address]: 1710, Building 11, Phase II of Tian'an Yungu Industrial Park, Gangtou Community, Bantian Street, Longgang District, Shenzhen,

[Postcode]: 518129 [Tel]: +86 - 755-89505533

[E-mail]: sales88@velpeau.com [Website]: www.velpeau.com

[Date of Preparation]: Jan.19,2017 [Modification Date]: Jul.15,2023

www.velpeau.com —